

ENJOYCARE

Sock Aid with Foam Handles

Model **EK-B105** · Joy Ours part no. **JOY-DALY-B105**


Put socks on while seated: stretch the sock over the form, step in and pull the foam-handled cords.

KEY FEATURES

- **Socks on while seated** Put your socks on sitting down, without reaching down to your feet.
- **Three-prong form** Holds the sock opening stretched so your foot can slide straight in.
- **Notched edges** Grip the sock so it stays on the form while you pull the cords.
- **Foam handles** Soft foam handles on the cord ends are easy to hold and pull.
- **Doubles as a shoehorn** Use the moulded body to help your heel slide into a shoe.
- **Light and wipeable** About 163 g, with a smooth plastic body that wipes clean.


SPECIFICATIONS

Width	9.7-12.3 cm
Colour	Blue
Weight	approx. 163 g
Material	PP, foam handles
Body Length	approx. 21 cm
Cord Length	approx. 75 cm
Overall Length	approx. 25.7 cm
Handle Length	approx. 10 cm
Side Height (Each End)	approx. 6.6 cm and 8.5 cm

WHO IT SUITS

- People who find bending forward or at the hip difficult
- Anyone who dresses sitting down and struggles to reach their feet
- People who also need help getting a heel into a shoe
- Carers encouraging someone to put on their own socks

GOOD TO KNOW

- Load the sock onto the three-prong end first so its opening is stretched, then hold both foam handles and place the form in front of your foot.
- Keep pulling the cords upward until the sock is fully on and the form comes out.
- The smooth PP body can be wiped clean.