

ENJOYCARE

Grace 2-in-1 Rollator & Transit Chair

Model **EKR936** · Joy Ours part no. **JOY-WALK-R936**


A rollator that converts to an attendant-pushed transit chair, with a footrest, reversible backrest and parking brakes.

KEY FEATURES

- **Walk or be pushed** Switch between a rollator and a transit chair pushed by a carer.
- **Reversible backrest** Reverses and folds, with four height positions.
- **Flip-up footrest** Supports the feet in chair mode and detaches for transport.
- **Parking brakes** Hand brakes with a parking function hold it still when stopped.
- **Seven handle heights** Overall height from 79 to 93.5 cm.
- **Compact fold** Folds to 63 × 24.5 × 46.5 cm, with a removable footrest.


SPECIFICATIONS

SWL	136 kg
Brakes	Hand-operated with Parking Function
Storage	Removable Under-seat Bag
Backrest	Reversible & Foldable
Footrest	Flip-up & Detachable
Net Weight	9.3 kg
Seat Width	41 cm / 16"
Rear Wheels	8"
Seat Height	50 cm / 19.5"
Folded Width	24.5 cm / 9.5"

Front Wheels	8"
Folded Height	46.5 cm / 18.5"
Folded Length	63 cm / 25"
Overall Width	60 cm / 23.5"
Frame Material	Aluminium Alloy
Overall Height	79-93.5 cm / 31-37"
Overall Length	102 cm / 40"
Handle Height Adjustment	7 Positions
Backrest Height Adjustment	4 Positions

WHO IT SUITS

- People who can walk some of the way but tire on longer outings
- Families or carers who want to push the user when walking stops
- Appointments, shopping centres and day trips
- Users up to 136 kg who suit a 41 cm wide seat

GOOD TO KNOW

- Apply the parking brakes before sitting down or switching between modes.
- At 9.3 kg it is heavier than the walking-only Enjoycare rollators (5.45-7.2 kg).
- In chair mode it is designed to be pushed by a carer.