

ENJOYCARE

Lightweight Aluminium Transit Wheelchair

Model **EW2000C** · Joy Ours part no. **JOY-WCHR-W20C**



A 10 kg aluminium transit wheelchair for a carer to push, with handle brakes and a 28 cm folded width.

KEY FEATURES

- **10 kg aluminium frame** Light for carers to lift, with an anodised, corrosion-resistant finish.
- **Folds to 28 cm** The dual cross-brace frame folds narrow for the car boot or a cupboard.
- **Attendant hand brakes** Carers control speed and stopping from the push handles; parking brakes hold the chair still.
- **Seat carry handles** Built-in handles on the seat make folding and lifting easier.
- **Breathable padded seat** 600D Oxford fabric with breathable mesh padding on the seat and backrest.
- **Puncture-proof tyres** Solid PVC front and PU rear tyres that do not go flat.

SPECIFICATIONS

SWL	100 kg
Net Weight	10 kg
Rear Tyres	Solid PU
Seat Depth	41 cm / 16"
Seat Width	43 cm / 17"
Upholstery	600D Oxford Fabric with Breathable Mesh Padding
Front Tyres	Solid PVC
Seat Height	44 cm / 17.5"
Folded Width	28 cm / 11"
Frame Finish	Anodised

Overall Width	63 cm / 25"
Braking System	Attendant Hand Brakes & Parking Brakes
Frame Material	6061F Aluminium Alloy
Overall Height	90 cm / 35.5"
Overall Length	96 cm / 38"
Backrest Height	42 cm / 16.5"
Rear Wheel Size	12"
Front Wheel Size	6"
Frame Tube Diameter	22 mm
Frame Tube Wall Thickness	2 mm

WHO IT SUITS

- People who are pushed by a carer for outings and appointments
- Carers who need a chair light enough to lift into a car boot
- Car trips and outings where storage space is tight
- Users up to 100 kg

GOOD TO KNOW

- This is a transit chair for a carer to push. For self-propelling, look at the EW2000H20 with 20-inch rear wheels.
- Put the parking brakes on before the person gets in or out.