

ENJOYCARE

Steerable Knee Walker with Shock Absorber

Model **EKR201P** · Joy Ours part no. **JOY-WALK-201P**


A steerable knee walker: rest your injured leg on the padded platform and push along with your other foot.

KEY FEATURES

- **Steerable T-bar handlebar** Turn the handlebar to steer around corners and through everyday spaces with more control.
- **Spring shock absorber** Helps soften bumps and vibration as you roll, with a contoured cushion supporting your lower leg.
- **Two height adjustments** Seven handlebar settings, and five knee platform settings from 46.5 to 56.5 cm.
- **Hand brake and bell** A hand-operated brake controls speed and stopping, and the bell sits on the handlebar.
- **Folding handlebar stem** Folds down to make the 9.6 kg walker easier to store and transport.
- **Bag and flip-up footrest** A removable bag carries your essentials, and the front footrest folds out of the way when not needed.

SPECIFICATIONS

SWL	100 kg
Net Weight	9.6 kg
Frame Material	Carbon Steel
Overall Height	80.5–101 cm / 31.5–40"
Rear Wheel Size	8"
Front Wheel Size	8"
Overall Dimensions	87.5 × 40.5 cm / 34.5 × 16"

Packing Dimensions	90 × 30 × 38 cm / 35.5 × 12 × 15"
Handlebar Adjustment	7 Positions
Knee Platform Height	46.5–56.5 cm / 18.5–22"
Knee Platform Adjustment	5 Positions
Knee Platform Dimensions	35.5 × 17 cm / 14 × 6.5"

WHO IT SUITS

- People recovering from a foot or ankle injury or surgery who need weight kept off that foot
- People who find crutches hard work on the hands and underarms
- Moving around the home and other suitable surfaces during recovery
- Users up to 100 kg who fit a knee platform height of 46.5–56.5 cm

GOOD TO KNOW

- Check the knee platform range (46.5–56.5 cm) and overall height (80.5–101 cm) against the user before ordering.
- The four 8-inch wheels are made for suitable surfaces, so think about the ground on your usual routes.
- Fold the handlebar stem down before lifting the walker into a car. It weighs 9.6 kg.